


Our Ready-to-Go Targets


trusted · reproducible · fast

WEE1

Labeled MST






LabelFree TSA (nanoDSF)




Screening highlights

- Binding and inhibition assays
- Up to 100 T_m s and 75 K_D s/day (method dependent)
- 1-2 weeks turnaround time
- K_m and IC_{50} determination *via* ADP-Glo Biochemical assays




Staurosporine



Off-the-shelf structures

- Down to 1.8 Å resolution
- 2-4 weeks turnaround time
- 24 structures delivered to date



- ✓ Protein production
- ✓ nanoDSF
- ✓ MST
- ✓ ADP-Glo

